HybriDetect 2T

Universal lateral flow dipstick for simultaneous detection of two different analytes (proteins, genomic amplicates) labelled with FITC, biotin and digoxigenine; development platform
English: 1-8

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

Materials Supplied, Storage and Stability

Component	Cat.-No.	Content	Preparation	Storage	Shelf Life
HybriDetect 2T Dipsticks: Membrane coated with biotin-ligand and polyclonal (goat) digoxigenine antibody; polyclonal (rabbit) anti-FITC antibody in gold conjugate	MGDS2	2 x 50 tests (blue lamination)	Ready-to use	2 –8°C; Containers must be always closed (moisture protection)!	until expiry date
HybriDetect 2T Assay Buffer: Tris-buffered saline	MGCB2	2 vials à 10 mL (yellow cap)	Ready-to use	2 –8°C	until expiry date

Material Safety Data Sheets are available on request (look as well www.milenia-biotec.de).

Materials Required

- Pipets
- Pipet tips (containing protective filters for PCR)
- Reaction tubes or 96-well microtiter plate

Necessary Development Work - Development Platform

Development of two solutions (A and B). Each solution must contain two different labeled detectors for the respective analyte.

Following conditions are necessary:

- 1) **Solution A:** Detectors must be labeled with:
 - FITC (fluorescein isothiocyanate)
 - Biotin
- 2) **Solution B:** Detectors must be labeled with:
 - FITC (fluorescein isothiocyanate)
 - Digoxigenine
- 3) Use about 100 µL fluid (sample material and analyte-specific solutions, A and B) for the assay procedure.
- 4) The provided buffers (MGCB2) may be used as a basic for these analyte-specific solutions.

Example of use:

20 µl sample material and 80 µl analyte-specific solutions with 5 min incubation time.

Notice:

Volumes, analyte-specific solutions, and incubation time are part of the individual test development.

A basic procedure for the detection of genomic amplicons is explained on page 7: "Assay performance PCR products." One amplification product should be biotinylated and the specific hybridization probe should be FITC labeled. The second amplification product should be labelled with digoxigenine and the specific hybridization probe should be FITC labeled. Various amplification procedures (Polymerase Chain Reaction or isothermal amplifications like LAMP or RPA) can be used.

Warnings and Precautions

- All reagents should be stored at 2 – 8° C in their original containers.
- Before use, bring all reagents to room temperature (18 – 28° C).
- The expiration date of all components must be observed.
- Protect dipsticks from humidity; Container must always be closed.
- Touch and label only the foil-covered areas of the dipsticks (labeling area).

The assay buffer contains an anti-microbial reagent; therefore avoid contact with skin and/or mucous membranes.

Method

Milenia[®] HybriDetect 2T is a ready-to-use, universal test strip (dipstick), which bases on the lateral flow technology using gold particles. The dipstick is designed to develop qualitative or quantitative rapid test systems for simultaneous detection of two different analytes such as proteins, antibodies, or gene amplifications. The user needs to develop two analyte-specific solutions, with following conditions. Solution A: Contains a first detector (e.g. antibody, antigen, specific probe) labeled with FITC and a second one (e.g. antibody, primer) labeled with biotin, and solution B: Contains a first detector (e.g. antibody, antigen, specific probe) labeled with FITC again and but a second one (e.g. antibody, primer) labeled with digoxigenine, (see “Common Testprinciple” on page 5).

The sample to be determined is mixed with the developed analyte-specific solutions, and then the dipstick is placed into this solution.


The complexed analyte (A)-labeled with FITC and biotin, binds first to the gold-labeled FITC-specific antibodies in the sample application area of the dipstick as well as the complexed analyte B (labeled with FITC and digoxigenine). The gold complexes A and B diffuse over the membrane by capillarity. Only the analyte captured gold particles will bound when they overflow the immobilized biotin-ligand molecules at the respective test band (test band A- analyt A, test band B- analyt B) and generate there a red-blue band over the time. Not-captured gold particles flow over the control band and will be fixed there by species-specific antibodies. With increasing incubation time, the formation of an intensely colored control band appears.

Additional available Products


Product Name	Order No.	Content	Description
Milenia [®] HybriDetect	MGHD 1	100 tests	Dipsticks with one test band
Milenia [®] HybriDetect HS	MGHT 1	40 tests	Test strips with one test band (in a housing)
Milenia [®] POCScan Reader	MSCAN 1	1 device	Interpretation device for HybriDetect products
Milenia [®] HybriDetect Software	MGSW	1 CD	Software for the documentation and/or quantification of results
Milenia [®] Carrier	MGCA 1	15 pieces	Holder for insertion of Milenia [®] Dipsticks in the Milenia [®] POCScan Reader

Common Testprinciple

1) Ready-to-use (analyte-independent) dipstick

2) Addition of sample, resp. analyte A and B-specific solutions


Symbols:

	FITC
	Biotin
	Digoxigenine
	Goldparticle
	anti-FITC Antibody
	Analyte A Detector 1 + 2
	Analyte A (Protein, antibody, genetic amplicon)
	Analyte B Detector 1 + 2
	Analyte B (Protein, antibody, genetic amplicon)
	Biotin-Ligand
	anti-Digoxigenine Antibody
	anti-rabbit Antibody

3) Incubation

4) Results:

Control Band Dipstick


In any case, the control band must be always visible!

It is a control function and can not be used to assess the test band. If the control band is not visible after the incubation period, the result is invalid! The test must be repeated with a new dipstick!


Interpretation of Results

There are three possibilities to interpret the Milenia®HybriDetect 2T dipstick results:

1) Qualitative by visual interpretation

2) Semi-quantitative, e.g. by interpretation cards


3) Quantitative interpretation with the Milenia® POCScan Reader and appropriate software

Assay Performance “PCR Products”

1. Take the required number of dipsticks out of the container and mark them.
2. For each sample to be analyzed pipet **100 µL** HybriDetect Assay Buffer or individual developed buffer into a reaction tube or a well of a microtiter plate.
3. Pipet **5 – 10 µl** of the hybridization product directly on the sample application area, or **alternatively**: add **5 – 10 µl** of the hybridization product into the solution of the reaction tube / well.
4. Place the dipsticks with the sample application area into the solution and incubate them e. g. for **5 - 15 minutes** in an upright position.
5. In the end of incubation period, remove dipsticks from assay solution and interpret test results immediately.


Notice:

- If a higher analytical sensitivity is required, it could be helpful to increase the volume of the PCR product.
- Volumes, analyte-specific solutions, and incubation time are always part of the individual test development.

Interpretation of “PCR” Results

Applying an internal PCR control:

While the biotin-labeled **amplicon A** is mainly used for the detection of a specific target sequence the digoxigenine-labeled **amplicon B** is provided as an internal PCR amplification control.

Test band A	Test band B	Control band	Interpretation
positive	positive	positive	<ul style="list-style-type: none"> control band is clearly visible, test run is valid amplicon A is detected (positive) internal PCR control is positive (amplicon B)
negative	positive	positive	<ul style="list-style-type: none"> control band is clearly visible, test run is valid amplicon A is not detected (negative) internal PCR control is positive (amplicon B)
positive	negative	negative	<ul style="list-style-type: none"> test runs were not valid because control band did not appear
negative	positive		
negative	negative		
negative	negative	positive	<ul style="list-style-type: none"> rapid test performed properly because control band is clearly visible PCR amplification did not work due (amplicon B is not detectable)
positive			

Trouble Shooting „PCR“

Problem	Possible cause(s)	Recommendation
Control band is not visible.	<ul style="list-style-type: none"> a) wrong or destroyed assay buffer b) expiration date of dipsticks is exceeded c) wrong storage conditions of dipsticks 	apply new (fresh) chemicals
Negative result with dipstick but clearly visible band in agarose gel	<ul style="list-style-type: none"> a) detection of an unspecific PCR product in agarose gel b) hybridization was not successful 	Check identity of PCR product by Southern blotting or sequence analysis Check conditions of hybridization.
Mineral oil	<ul style="list-style-type: none"> a) mineral oil affects flow characteristics of the assay b) development of band might be hampered 	Remove PCR product very slowly from the bottom of the reaction vial.

Assay Sensitivity “PCR”

Analytical sensitivity of the Milenia® HybriDetect 2T is equivalent to agarose gel electrophoresis and subsequent staining with Ethidiumbromide. Independent on the size of the PCR product and the number of amplification cycles as low as 5 pg DNA could be detected.

Literature-References

Kiatpathomchai W, et al, J Virol Methods (2008); 153: 214-217
 Puthawibool T, et al, J Virol Methods (2009); 156 (1-2): 27-31
 Jaroenram W, et al, Mol Cell Probes (2009); 23 (2): 65-70
 Nimitphak T, et al, Mol Cell Probes (2009); 24 (1): 1-5
 Kikuchi T, et al, Nematology (2009) ; 99 (12): 1365-1369
 Piepenburg O, et al, PloS Biology 2006, Volume 4, Issue 7, e204

HybriDetect 2T

Universal lateral flow dipstick for simultaneous detection of two different analytes (proteins, genomic amplicates) labelled with FITC, biotin and digoxigenine; development platform

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